


MIS II


Posterior Spinal Fixation System for Minimally Invasive Surgery

“Simplified and versatile fixation system for MIS”



- **Long Arm Design**

- Designed for percutaneous minimally invasive surgery, minimize trauma to muscle and soft tissue
- Provide 20 mm of controlled rod reduction



- **Wider Angle**

- Offer 40° of angulation in any direction for intraoperative flexibility (L5-S1 is available)



- **Progressive Double Thread**

- Enhance bone purchase and pull-out resistance
- Promote rapid screw placement



Cannulated Screw



Bone Cement Screw


- **Curve Smooth Rod**




Bullet-nosed design to ease navigation through soft tissue

- **Rod Holder for MIS Surgery**

- Simplified and streamlined rod placement procedure
- Provide more flexible percutaneous rod rotation control



Elpis Multi MISS Screw



Cannulated Screw

REF	Diameter (mm)	Length (mm)	REF	Diameter (mm)	Length (mm)
E1186-4525	4.5	25	E1186-6525	6.5	25
E1186-4530	4.5	30	E1186-6530	6.5	30
E1186-4535	4.5	35	E1186-6535	6.5	35
E1186-4540	4.5	40	E1186-6540	6.5	40
E1186-4545	4.5	45	E1186-6545	6.5	45
E1186-4550	4.5	50	E1186-6550	6.5	50
E1186-4555	4.5	55	E1186-6555	6.5	55
E1186-4560	4.5	60	E1186-6560	6.5	60
E1186-4565	4.5	65	E1186-6565	6.5	65
E1186-4570	4.5	70	E1186-6570	6.5	70
E1186-4575	4.5	75	E1186-6575	6.5	75
E1186-4580	4.5	80	E1186-6580	6.5	80
E1186-5025	5.0	25	E1186-7025	7.0	25
E1186-5030	5.0	30	E1186-7030	7.0	30
E1186-5035	5.0	35	E1186-7035	7.0	35
E1186-5040	5.0	40	E1186-7040	7.0	40
E1186-5045	5.0	45	E1186-7045	7.0	45
E1186-5050	5.0	50	E1186-7050	7.0	50
E1186-5055	5.0	55	E1186-7055	7.0	55
E1186-5060	5.0	60	E1186-7060	7.0	60
E1186-5065	5.0	65	E1186-7065	7.0	65
E1186-5070	5.0	70	E1186-7070	7.0	70
E1186-5075	5.0	75	E1186-7075	7.0	75
E1186-5080	5.0	80	E1186-7080	7.0	80
E1186-5525	5.5	25	E1186-7525	7.5	25
E1186-5530	5.5	30	E1186-7530	7.5	30
E1186-5535	5.5	35	E1186-7535	7.5	35
E1186-5540	5.5	40	E1186-7540	7.5	40
E1186-5545	5.5	45	E1186-7545	7.5	45
E1186-5550	5.5	50	E1186-7550	7.5	50
E1186-5555	5.5	55	E1186-7555	7.5	55
E1186-5560	5.5	60	E1186-7560	7.5	60
E1186-5565	5.5	65	E1186-7565	7.5	65
E1186-5570	5.5	70	E1186-7570	7.5	70
E1186-5575	5.5	75	E1186-7575	7.5	75
E1186-5580	5.5	80	E1186-7580	7.5	80
E1186-6025	6.0	25	E1186-7585	7.5	85
E1186-6030	6.0	30	E1186-7590	7.5	90
E1186-6035	6.0	35	E1186-7595	7.5	95
E1186-6040	6.0	40	E1186-75100	7.5	100
E1186-6045	6.0	45	E1186-75105	7.5	105
E1186-6050	6.0	50	E1186-75110	7.5	110
E1186-6055	6.0	55	E1186-8025	8.0	25
E1186-6060	6.0	60	E1186-8030	8.0	30
E1186-6065	6.0	65	E1186-8035	8.0	35
E1186-6070	6.0	70	E1186-8040	8.0	40
E1186-6075	6.0	75	E1186-8045	8.0	45
E1186-6080	6.0	80	E1186-8050	8.0	50
E1186-6225	6.25	25	E1186-8055	8.0	55
E1186-6230	6.25	30	E1186-8060	8.0	60
E1186-6235	6.25	35	E1186-8065	8.0	65
E1186-6240	6.25	40	E1186-8070	8.0	70
E1186-6245	6.25	45	E1186-8075	8.0	75
E1186-6250	6.25	50	E1186-8080	8.0	80
E1186-6255	6.25	55	E1186-8085	8.0	85
E1186-6260	6.25	60	E1186-8090	8.0	90
E1186-6265	6.25	65	E1186-8095	8.0	95
E1186-6270	6.25	70	E1186-80100	8.0	100
E1186-6275	6.25	75	E1186-80105	8.0	105
E1186-6280	6.25	80	E1186-80110	8.0	110

Multiaxial MISS Bone Cement Screw



Bone Cement Screw



Pusher for
Cement Delivery

REF	Diameter (mm)	Length (mm)	REF	Diameter (mm)	Length (mm)
E1326-4525M	4.5	25	E1326-6525M	6.5	25
E1326-4530M	4.5	30	E1326-6530M	6.5	30
E1326-4535M	4.5	35	E1326-6535M	6.5	35
E1326-4540M	4.5	40	E1326-6540M	6.5	40
E1326-4545M	4.5	45	E1326-6545M	6.5	45
E1326-4550M	4.5	50	E1326-6550M	6.5	50
E1326-4555M	4.5	55	E1326-6555M	6.5	55
E1326-4560M	4.5	60	E1326-6560M	6.5	60
E1326-4565M	4.5	65	E1326-6565M	6.5	65
E1326-4570M	4.5	70	E1326-6570M	6.5	70
E1326-4575M	4.5	75	E1326-6575M	6.5	75
E1326-4580M	4.5	80	E1326-6580M	6.5	80
E1326-5025M	5.0	25	E1326-7025M	7.0	25
E1326-5030M	5.0	30	E1326-7030M	7.0	30
E1326-5035M	5.0	35	E1326-7035M	7.0	35
E1326-5040M	5.0	40	E1326-7040M	7.0	40
E1326-5045M	5.0	45	E1326-7045M	7.0	45
E1326-5050M	5.0	50	E1326-7050M	7.0	50
E1326-5055M	5.0	55	E1326-7055M	7.0	55
E1326-5060M	5.0	60	E1326-7060M	7.0	60
E1326-5065M	5.0	65	E1326-7065M	7.0	65
E1326-5070M	5.0	70	E1326-7070M	7.0	70
E1326-5075M	5.0	75	E1326-7075M	7.0	75
E1326-5080M	5.0	80	E1326-7080M	7.0	80
E1326-5525M	5.5	25	E1326-7525M	7.5	25
E1326-5530M	5.5	30	E1326-7530M	7.5	30
E1326-5535M	5.5	35	E1326-7535M	7.5	35
E1326-5540M	5.5	40	E1326-7540M	7.5	40
E1326-5545M	5.5	45	E1326-7545M	7.5	45
E1326-5550M	5.5	50	E1326-7550M	7.5	50
E1326-5555M	5.5	55	E1326-7555M	7.5	55
E1326-5560M	5.5	60	E1326-7560M	7.5	60
E1326-5565M	5.5	65	E1326-7565M	7.5	65
E1326-5570M	5.5	70	E1326-7570M	7.5	70
E1326-5575M	5.5	75	E1326-7575M	7.5	75
E1326-5580M	5.5	80	E1326-7580M	7.5	80
E1326-6025M	6.0	25	E1326-7585M	7.5	85
E1326-6030M	6.0	30	E1326-7590M	7.5	90
E1326-6035M	6.0	35	E1326-7595M	7.5	95
E1326-6040M	6.0	40	E1326-75100M	7.5	100
E1326-6045M	6.0	45	E1326-75105M	7.5	105
E1326-6050M	6.0	50	E1326-75110M	7.5	110
E1326-6055M	6.0	55	E1326-8025M	8.0	25
E1326-6060M	6.0	60	E1326-8030M	8.0	30
E1326-6065M	6.0	65	E1326-8035M	8.0	35
E1326-6070M	6.0	70	E1326-8040M	8.0	40
E1326-6075M	6.0	75	E1326-8045M	8.0	45
E1326-6080M	6.0	80	E1326-8050M	8.0	50
E1326-6225M	6.25	25	E1326-8055M	8.0	55
E1326-6230M	6.25	30	E1326-8060M	8.0	60
E1326-6235M	6.25	35	E1326-8065M	8.0	65
E1326-6240M	6.25	40	E1326-8070M	8.0	70
E1326-6245M	6.25	45	E1326-8075M	8.0	75
E1326-6250M	6.25	50	E1326-8080M	8.0	80
E1326-6255M	6.25	55	E1326-8085M	8.0	85
E1326-6260M	6.25	60	E1326-8090M	8.0	90
E1326-6265M	6.25	65	E1326-8095M	8.0	95
E1326-6270M	6.25	70	E1326-80100M	8.0	100
E1326-6275M	6.25	75	E1326-80105M	8.0	105
E1326-6280M	6.25	80	E1326-80110M	8.0	110

Curve Smooth Rod



REF	Diameter (mm)	Length (mm)
R1600-5504C	5.5	40
R1600-55045C	5.5	45
R1600-5505C	5.5	50
R1600-55055C	5.5	55
R1600-5506C	5.5	60
R1600-55065C	5.5	65
R1600-5507C	5.5	70
R1600-55075C	5.5	75
R1600-5508C	5.5	80
R1600-55085C	5.5	85
R1600-5509C	5.5	90
R1600-55095C	5.5	95
R1600-5510C	5.5	100
R1600-5511C	5.5	110
R1600-5512C	5.5	120
R1600-5513C	5.5	130
R1600-5514C	5.5	140
R1600-5515C	5.5	150
R1600-5516C	5.5	160
R1600-5517C	5.5	170
R1600-5518C	5.5	180
R1600-5519C	5.5	190
R1600-5520C	5.5	200

Inside Hex Set Screw



REF	Height (mm)
A1110-1060	5.0

Pusher



REF
PD-T101-A